


# Monitect®

## Drug Screen Cassette Test



### Drugs of Abuse Test Cut-off Levels:

Amphetamine (AMP)	1000 ng/mL
Barbiturates (BAR)	300 ng/mL
Benzodiazepines (BZO)	300 ng/mL
Cocaine (COC)	300 ng/mL
Ecstasy (MDMA)	500 ng/mL
Marijuana (THC)	50 ng/mL
Methamphetamine (MET)	1000 ng/mL
Methadone (MTD)	300 ng/mL
Opiates (OPI)	300/2000 ng/mL
Oxycodone (OXY)	100 ng/mL
Phencyclidine (PCP)	25 ng/mL
Tricyclic Antidepressants (TCA)	1000 ng/mL

### New:

Buprenorphine* (BUP)	5 ng/mL
----------------------	---------

\*For Sale Outside US Only

### Tests for Adulterants:

Oxidizing Agents  
pH  
Creatinine  
Nitrite

### Coming Soon:

Specific Gravity

- Simultaneous detection of multiple drugs and adulterants
- Available with or without adulteration tests
- Single or multiple test configuration
- User-friendly, reliable, accurate
- One-step qualitative assay
- Easy interpretation of results (5 minutes or less)
- FDA Cleared



